


MRS-6 Retainer


July 2018



Two views of the MRS-6 retainer (SS304) are shown above.

The top section holds the MRS-6 and is 12.5 mm X 3 mm. The bottom is a 3 mm aluminum pin. The MRS-6 is approx. 0.2 mm below the top surface secured by a "C" clip from below.

Why purchase the retainer? The MRS-6 is a thin 2 X 2 mm silicon chip which is fragile. It should be handled with a vacuum or plastic tweezer to prevent edge chipping and possible cleaving.



GELLER
MICROANALYTICAL
LABORATORY

426e BOSTON ST., TOPSFIELD, MA 01983-1216
TEL 978 887-7000 FAX 978-887-6671
<http://www.gellermicro.com>